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en

Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION

Chinese rubber cups are an ideal tool for home vacuum massage. Thanks to their different sizes, they can be adjusted to different areas of the body, such as the face, arms, thighs or abdomen. The cups are easy to use, flexible and durable, making them an excellent choice for people who want to take care of their body at home.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

INSTRUCTIONS FOR USE

- **Preparation:**
 - Before first use, rinse the brush under running water.
 - Make sure the surface you want to attach the brush to is clean and dry.
- **Installation:**
 - Place the brush on the surface of the bathtub or shower.
 - Press firmly to make sure the suction cups adhere properly.
- **Use :**
 - Dampen the brush with water.
 - Apply a little soap or bath gel to your hair.
 - Gently move your foot over the brush in a back and forth motion.
 - The brush automatically massages and cleanses the skin of the feet.
- **After use:**
 - Rinse the brush under running water to remove any soap scum and dirt.
 - Detach the brush by pulling on one of the corners.
 - Leave to dry in a dry place.

TECHNICAL DATA

- Material: High quality rubber.
- Cup sizes: 25mm, 40mm, 50mm, 65mm.
- Color: Blue.
- Application: Vacuum massage, improved circulation, cellulite reduction.

SAFETY INSTRUCTIONS

- Before installing and using the product, please read this instruction manual carefully.
- This product may be used by children aged 12 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Children must not play with the product.
- Do not use the brush on damaged or irritated skin.
- Do not use in case of problems with varicose veins, broken capillaries or fragile skin.
- Cleaning and maintenance may not be performed by children without supervision.

- Massage on one body part should not exceed 10 minutes.
- Wash cups after each use to avoid bacteria build-up.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Keep the foil packaging of the product away from children (risk of suffocation)
- Avoid sources of ignition: avoid contact with flames or heat-generating devices (may cause fire)
- Protect the product from extreme temperatures, direct sunlight, high humidity, flammable gases, solvents.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is prohibited to use the product if any part is damaged.
- Always use the product as intended.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

The product meets European and national requirements on the safety of devices and products.
We reserve the right to make changes to text, design and product data without notice.

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